

LIVE WEBINAR

# Unlocking the Inner Ecosystem

Decoding the Gut Microbiome for Clinical Application

Wednesday, Mar 29 | 7 - 8PM AEDT

**CO-BIOME**  
by MICROBA

1

## Meet your speakers



**Hayley Parcell**  
Nutritionist and Head of Co-Biome™ Healthcare



**Dr Brad Leech**  
Nutritionist and Lead Clinical Educator



**Krystyna Sullivan**  
Naturopath and Clinical Application Specialist



All participants have been muted



There is an optional 15 minutes for questions at the end



Add your questions in the chat and we will come back to them at the end

**CO-BIOME**  
by MICROBA

2

## Acknowledgement of Country

## Disclaimers

- The information provided in this webinar is for the use of qualified healthcare professionals.
- The information contained in this webinar is in no way to be taken as prescriptive or to replace a healthcare professional's duty of care and personalised care practices.
- The clinical opinions and patient case studies shared by presenters are solely those of the individual presenters and do not necessarily represent the view of Co-Biome.

## What we'll be covering

---

Introducing Co-Biome™

---

The role of the gut microbiome in clinical practice

---

Decoding the gut microbiome for clinical application

---

The Co-Biome™ MetaXplore™ range

---

Case study: Interpreting a MetaXplore™ GI Plus report

---

Q & A

CO-BIOME  
by MICROBA

5

## Introducing Co-Biome™

Hayley Parcell

CO-BIOME  
by MICROBA



6



## Improving human health through precision microbiome testing

Co-Biome™ believes that the future of good health lies within us and through accurately unlocking the complexity of the gut microbiome, we can better manage patient health for the future.

7

### What's the Co-Biome difference?

<b>Contribution</b> We are contributing to improve human health	<b>Connection</b> We are connecting with healthcare professionals and their patients	<b>Commitment</b> We are committed to clinically useful information	<b>Courage</b> We show courage to pioneer evidence-based solutions	<b>Conscious</b> We are conscious of our impact
Whole gut microbiome analysis through metagenomic assessment of over 28,000 species	Supportive team of clinical application specialists	Developed by a team of clinicians, specifically for clinicians	Scientifically graded statements on the evidence	Local Australian NATA accredited lab with world-class technology




8



# The role of the gut microbiome in clinical practice

Dr Brad Leech


CO-BIOME  
MICROBIA



9

## What is the Microbiome


The microbiome is defined as “a characteristic microbial community occupying a reasonably well-defined habitat which has distinct physio-chemical properties” (*Whipps et al, 1988*)



CO-BIOME  
MICROBIA

10

**The colonic mucus provides a home and a fuel source for the microbiome**




CO-BIOME  
MICROBIA

11

**The microbes within the gut microbiome consume:**

- dietary fibre and resistant starch
- excess dietary protein
- host proteins
- mucin glycans




CO-BIOME  
MICROBIA


12

### The amount of metabolite produced by the microbiome reflects:


- the number of microbial genes that produce different metabolites
- availability of fuel sources
- environmental conditions in the gut



### How can the gut microbiome influence health?



## What can influence a patient's microbiome health?



CO-BIOME  
by MICROBA

15

## Decoding the gut microbiome for clinical application

Dr Brad Leech


CO-BIOME  
by MICROBA



16



## gastrointestinal health



- Gastrointestinal disorders, including irritable bowel syndrome (IBS), intestinal permeability, and inflammatory bowel disease (IBD)
- Digestive complaints, including constipation, diarrhoea, bloating and abdominal pain
- Hormonal imbalance
- Metabolic and weight concerns
- Immune system health concerns
- Chronic inflammation

CO-BIOME  
by MICROBA

17

## How can gut microbiome analysis help us?


Microbiome's contribution to health	Dietary direction	Supplements	Patient compliance
<p>How is my patient's microbiome contributing to their clinical presentation?</p> <p>What is the severity of my patient's gut health?</p> <p>How is my patient's microbiome contributing to health or disease risk?</p>	<p>What personalised dietary advice should I provide my patient?</p> <p>How is my patient's diet impacting their microbiome?</p>	<p>What are the prebiotic, probiotic and polyphenol supplements which should be used or avoided for my patient?</p> <p>Which supplements should I prioritise?</p>	<p>Allow us to monitor and support patient compliance and treatment over time.</p>

CO-BIOME  
by MICROBA

**"Test don't guess"**

18

## Decoding the gut microbiome for clinical application



CO-BIOME  
MICROBIA


19

## 6 key questions to answer

1. What is the microbial diversity and richness?

2. What is the whole microbiome's capacity to produce or consume metabolites associated with health?

6. Are there alterations in gut function or environment?



3. Are pathogenic bacteria or protist parasites present? (if required)

4. Who makes up the whole ecosystem – bacteria, archaea, fungi, protists?


5. Which species are present and what is their abundance?

CO-BIOME  
MICROBIA


20

## 1. What is the microbial diversity and richness?


### DIVERSITY



**LOW DIVERSITY:**  
Low richness  
High evenness



**LOW DIVERSITY:**  
High richness  
Low evenness



**HIGH DIVERSITY:**  
High richness  
High evenness

A healthy person's sample will contain between 110 – 244 species.<sup>1</sup>

Higher microbial diversity is predictive of a more stable microbiome.<sup>2</sup>

Low diversity may be associated with increased systemic inflammation.<sup>3, 4</sup>


CO-BIOME  
MICROBA

1. Data from Microba's Future Insights program. 2. Byrd, 2021 3. Mokkala, 2020, 4. Zhernakova, 2016

21

## 2. What is the microbiome's capacity to produce or consume metabolites associated with health?

### FUNCTION



CO-BIOME  
MICROBA

22

## 2. What is microbiome's capacity to produce or consume metabolites associated with health?

FUNCTION

	MICROBIAL MARKER	CLINICAL RELEVANCE
Short chain fatty acids	Butyrate	Main fuel source for gut cells. Low levels associated with intestinal and systemic inflammation and impaired intestinal barrier integrity.
	Acetate	Most abundant SCFA produced in the gut. Can be converted by some species to butyrate.
	Propionate	Low levels associated with intestinal inflammation. SCFA produced by the gut microbiome.
		Optimal production associated with normal gut transit time and immune balance within the GIT.
Health indicators	3-indolepropionic acid (IPA)	Beneficial substance produced by some gut bacteria when they break down tryptophan. Low levels may be associated with intestinal and systemic inflammation and impaired barrier integrity.
	Mucin	Component of the mucous layer that lines and protects the epithelial gut barrier. High levels of mucin consuming microbes may be associated with intestinal inflammation.
	Oxalate consumption	Key component of calcium oxalate kidney stones. Low microbial oxalate consumption associated with increased urinary oxalate excretion. May be reduced in patients with recurrent kidney stones.
	Hydrogen sulphide	Produced by gut microbes when they break down sulphur-containing compounds. It is responsible for the rotten egg smell of flatulence. Can be protective of the gut at low levels. At high levels can disrupt intestinal barrier integrity.
	Methane	Gas produced by some species of the gut microbiome, primarily through reduction of carbon dioxide and hydrogen. Elevated levels associated with increased intestinal transit time and constipation.
	<i>B. fragilis</i> toxin	Normal inhabitant of the human gut. Small proportion can secrete a toxin called fragilysin. High levels may impair intestinal barrier integrity.
	Beta-glucuronidase	Bacterial enzyme that can re-activate a wide variety of drugs and hormones. High levels may affect drug response and toxicity.
	Branched-chain amino acids (BCAAs)	Derived from the diet and produced by the gut microbiome which can contribute to elevated plasma BCAAs levels. High levels associated with systemic inflammation. Muscle important in regulating BCAA levels through regular physical activity.
	Hexa-acylated lipopolysaccharide (hexa-LPS)	Inflammatory compound produced by some species of bacteria within the Proteobacteria phylum. High levels associated with intestinal and systemic inflammation and impaired intestinal barrier integrity, via activation of immune receptor TLR4.
	Trimethylamine (TMA)	Compound produced by some microbes from breakdown of choline and carnitine. TMA produced in gut transported to liver, where it can convert to TMAO. High level of TMAO in blood plasma is strongly associated with cardiometabolic disease.

CO-BIOME  
+MICROBIA

Green = positively associated

Grey = scale dependent

Red = negatively associated

23

## Clinical Application of Hexa-LPS

### Clinical Relevance

Hexa-acylated lipopolysaccharides (hexa-LPS) are bacterial cell wall components of bacteria within the Gammaproteobacteria class.

High hexa-LPS may be associated with intestinal and systemic inflammation and impaired intestinal barrier integrity.



### Treatment Considerations

Clinicians may want to prescribe galacto-oligosaccharides (GOS) or the probiotic combination of *Lactobacillus gasseri* KS-Y3, *Bifidobacterium bifidum* G9-Y and *Bifidobacterium longum* MM2 to reduce hexa-LPS.

An increased omega-3 to saturated fat ratio in the diet may reduce blood levels of hexa-LPS after meals.

CO-BIOME  
+MICROBIA

Lyte, 2016; Ahola, 2017

24



## Clinical Application of IPA

### Clinical Relevance

3-indolepropionic acid (IPA) is a beneficial substance produced by some gut bacteria when they break down the amino acid tryptophan.

Low levels of IPA may be associated with intestinal and systemic inflammation and impaired intestinal barrier integrity.



### Treatment Considerations

Clinicians may want to suggest foods rich in ellagic acid (e.g., chestnuts, blackberries, and ellagic acid enriched pomegranate juice), as well as a Mediterranean diet for patients with low IPA production.

25

### 3. Are pathogenic bacteria or protist parasites present?

#### Clinical consequence

Nausea, vomiting, diarrhoea, abdominal discomfort, bloody stool, gas.<sup>1</sup>

Impacts gut environment.<sup>2</sup>

#### Exposure


Contaminated food, drinking water, animals, waste, soil.

#### Pathogenic bacteria

*Yersinia enterocolitica*  
*Aeromonas* spp.  
*Campylobacter* spp.  
*Salmonella* spp.  
*Vibrio* spp.  
*Clostridium difficile* toxin B  
Hypervirulent *Clostridium difficile*  
*E. coli* pathotypes x 6

#### Protist parasites


*Giardia lamblia*  
*Entamoeba histolytica* *Cryptosporidium* spp.  
*Dientamoeba fragilis*  
*Cyclospora cayatanensis*



26

#### 4. Who makes up the whole ecosystem – bacteria, archaea, fungi, protists?

SPECIES



Shotgun metagenomics can report on all species that make up the whole gut microbiome community<sup>^</sup>

It's important to consider individual species in the context of the whole microbiome


CO-BIOME  
MICROBA

<sup>^</sup> Shotgun metagenomics can report on all species with a relative abundance above 0.01% including non-diagnostic eukaryotes (fungi and protist parasites).

27

#### 5. What species are present and what is their abundance?

SPECIES



Data from Microba case study.

There are no species a healthy microbiome must contain

The most common species, *Blautia\_A wexlerae*, is found in 98% of healthy cohort<sup>1</sup>

The amount of microbes doing the job is more important than which species are doing the job

CO-BIOME  
MICROBA

1. Data from Microba healthy cohort.

28


## SPECIES



15

## 6. Are there alterations of gut function or environment?

GUT FUNCTION & ENVIRONMENT



CO-BIOME  
-MICROBA

31

## Clinical Application of Functional Markers

GUT FUNCTION & ENVIRONMENT

Functional Marker	Clinical indication	Clinical relevance
Calprotectin	• Acute intestinal inflammation and estimated degree of inflammation	• Distinguish active IBD from IBS • Monitor disease activity and relapse in IBD and colorectal cancer
Lactoferrin	• Intestinal inflammatory activity	• Monitor disease activity /treatment /relapse in IBD
Faecal pH	• Gut transit time	• Elevated pH may indicate longer transit time • Low pH may indicate rapid transit time • Seen in patients with lactose intolerance • Estimates absorption of short-chain fatty acids (SCFAs)
Faecal Occult Blood	• Intestinal bleeding	• Early diagnosis can significantly reduce risk of a serious colorectal disease
Pancreatic elastase	• Pancreatic exocrine function	• Diagnosis or exclusion of exocrine pancreatic insufficiency • Monitor exocrine pancreatic function in cystic fibrosis, diabetes mellitus, chronic pancreatitis
Zonulin	• Increased intestinal permeability	• Allows substances from gut lumen to pass across the epithelium and activate immune reactions
Secretory IgA	• Intestinal inflammation and increased intestinal permeability	• Major role in preventing adherence of microbes to mucosal sites, activation of the alternative complement pathway and inflammatory reactions

CO-BIOME  
-MICROBA

32



## Key Take Aways


1. Assess microbial diversity and richness

2. Assess the whole microbiome's capacity to produce or consume metabolites associated with health

*butyrate, acetate, propionate, 3-indolepropionic acid (IPA) hexa-acylated lipopolysaccharide (hexa-LPS), trimethylamine (TMA), hydrogen sulphide, branched-chain amino acids (BCAA), B.fragilis toxin, methane, beta-glucuronidase*

6. Assess alterations of gut function or environment

*faecal calprotectin, faecal occult blood, faecal pH, lactoferrin, pancreatic elastase, secretory IgA, zonulin*



3. Identify pathogenic bacteria and protist parasite presence using qPCR when necessary

4. Assess the WHOLE microbiome using metagenomics

5. Understand the presence and abundance at a species-level view using metagenomics


CO-BIOME  
MICROBIA

33

## The Co-Biome MetaXplore range

Krystyna Sullivan

CO-BIOME  
MICROBIA



34

## Introducing the MetaXplore™ range



### MetaXplore™

MetaXplore™ provides a metagenomic driven gut microbiome profile, together with the latest research insights for healthcare professionals.  
Technology: metagenomics

\$369



### MetaXplore™ GI

MetaXplore™ GI provides the same comprehensive microbiome profile as MetaXplore™ as well as reporting on seven gastrointestinal health markers and science backed clinical insights to assist clinical decision-making and intervention.  
Technology: metagenomics + diagnostic GI health markers + faecal pH

\$489



### MetaXplore™ GI Plus

MetaXplore™ GI Plus is Co-Biome's most comprehensive functional gut microbiome profile. It provides all the features found in MetaXplore™ and MetaXplore™ GI, plus targeted bacteria and parasite panels using quantitative polymerase chain reaction (qPCR) technology for diagnostic pathogen detection.  
Technology: metagenomics + diagnostic GI health markers + faecal pH + qPCR

\$529

CO-BIOME  
+MICROBA

35

### MetaXplore™ range test indications:

- Gastrointestinal disorders, including irritable bowel syndrome (IBS), intestinal permeability and inflammatory bowel disease (IBD)
- Digestive complaints, including constipation, diarrhoea, bloating and abdominal pain
- Hormonal imbalance
- Metabolic and weight concerns
- Immune system health concerns
- Chronic inflammation

### What does the MetaXplore™ range measure?

- Diversity
- Species
- Function
- Gut function & environment<sup>1</sup>



### What insights can the MetaXplore™ range provide?

The MetaXplore™ range can show how the gut microbiome and gastrointestinal health is impacting health categories including:



Intestinal inflammation



Systemic inflammation



Detox / retox



Intestinal motility



Digestive secretions









Intestinal barrier

CO-BIOME  
+MICROBA

36

## How is their ecosystem functioning?

Gut microbiome and gastrointestinal health impact on health categories

	 Intestinal inflammation	 Digestive secretions	 Intestinal motility	 Systemic inflammation	 Detox / retox	 Intestinal barrier
MICROBIOME FUNCTION	Acetate Butyrate Hexa-LPS IPA Mucin consumption Propionate		Diversity Methane Propionate Richness	BCAAs Butyrate Diversity Hexa-LPS IPA Richness Trimethylamine	Beta-glucuronidase Oxalate consumption	<i>B.fragilis</i> toxin Hexa-LPS IPA Hydrogen sulphide
GUT FUNCTION & ENVIRONMENT (optional)	Calprotectin Lactoferrin Occult blood Secretory IgA	Pancreatic elastase	Faecal pH			Secretory IgA Zonulin

CO-BIOME  
+MICROBA




37

The MetaXplore™ range provides **personalised clinical and research insights** using scientifically graded statements on the evidence for **diet** and **lifestyle** interventions **plus probiotic, prebiotic, nutrient** and **polyphenol** supplementation

CLINICAL INSIGHT	Grades / Code	Description
To reduce faecal calprotectin, advising patients to follow a Mediterranean diet may be considered.	A	Body of evidence can be trusted to guide practice
EVIDENCE D SHOW REFERENCES	B	Body of evidence can be trusted to guide practice in most situations
	C	Body of evidence provides some support for recommendation, but care should be taken in its application
RESEARCH INSIGHT	D	Body of evidence is weak, and recommendation must be applied with caution
Carnitine supplementation has been shown to increase plasma trimethylamine N-oxide (TMAO). When aiming to reduce plasma TMAO, limiting or avoiding carnitine supplementation may be effective.	PP, H	Body of evidence is observational only and must be applied with caution
EVIDENCE A SHOW REFERENCES	PP, IV	Body of evidence is in vitro and must be applied with a high degree of caution

CO-BIOME  
+MICROBA

38

TEST	WHEN TO REFER
<b>MetaXplore™</b> 	<ul style="list-style-type: none"> <li>Easiest swab collection and return process</li> <li>Reason for test is to explore role of microbiome in: <ul style="list-style-type: none"> <li>systemic inflammation, e.g. cardiovascular disease, obesity</li> <li>detox/ retox, e.g. hormone imbalance, oxalate processing</li> </ul> </li> </ul>
<b>MetaXplore™ GI</b> 	<ul style="list-style-type: none"> <li>Provides combination of microbial markers to assess impact of microbiome, and gastrointestinal health markers to assess gut function and environment</li> <li>Reason for test may be related to: <ul style="list-style-type: none"> <li>altered intestinal motility, e.g. constipation, diverticular disease</li> <li>compromised intestinal barrier function, e.g. autoimmune disease, food sensitivities</li> <li>intestinal inflammation, e.g. Inflammatory bowel disease, allergy</li> <li>pancreatic function, e.g. weight loss, alcoholism</li> </ul> </li> </ul>
<b>MetaXplore™ GI Plus</b> 	<ul style="list-style-type: none"> <li>Most comprehensive test available</li> <li>Provides qPCR diagnostic panel for pathogenic bacteria and protist parasites</li> <li>Recommended for patients: <ul style="list-style-type: none"> <li>wanting comprehensive gold standard test</li> <li>with symptoms of pathogen infection, e.g. abdominal pain, blood in stool, travel history</li> <li>who may have post-infective irritable bowel syndrome</li> </ul> </li> </ul>


CO-BIOME  
by MICROBA

39

## Case study: Interpreting MetaXplore™ GI Plus Reports

Dr Brad Leech

CO-BIOME  
by MICROBA



40



**Age** - 52

**BMI** - 27.5 (overweight)

**Diagnosis** - Hashimoto's disease, systemic lupus erythematosus, Sjögren's syndrome, Raynaud's syndrome, fibromyalgia, migraine headache

**Gastrointestinal** - nausea, food sensitivities, burning indigestion pain, chronic diarrhoea, reflux, intermittent constipation

**Family history** - heart disease, colon cancer

**Systemic** - chronic pain, high blood pressure

**Dietary** - low carbohydrate diet, low fibre intake

**Medication** - Nexium 40mg, Gaviscon b.d.

## Patient Case Study

Chronic autoimmune disease and gastrointestinal symptoms



41

## Patient Centred Care

### Important considerations when using research and clinical insights

- The patient's clinical presentation including:
  - Symptoms
  - Medical history
  - Physical examination
- Additional pathology results
- Any patient sensitivities or preferences
- Important areas of the report to address first
- Any conflicting insights suggested

Based on the patient's clinical presentation and MetaXplore™ results the clinician should determine the most appropriate course of treatment.

48

## Patient treatment plan

Supplement	Dosage	Duration	Related condition
Fish oil with SPM	Two with breakfast and lunch	3 months	Pain and inflammation
Curcumin	Two with breakfast and dinner	2 months	Intestinal and systemic inflammation
PEA	One with breakfast and dinner	2 months	Pain and inflammation

Dietary/ Lifestyle	Related condition
High plant-based diet	Low diversity, high TMAO
Breakfast: rolled oats with blackberries and walnuts	Low butyrate, low IPA (contains beta glucan) and ellagic acid
Increase fermented foods	Low microbial diversity
Weekly gardening	Low microbial diversity
Green banana flour	Low butyrate




49

## Key takeaways

- Assessing a patient's gut microbiome and gastrointestinal health should be a consideration for many of your patients
- Consider diversity, function, species and gut function and environment when assessing a patient's gut microbiome and gastrointestinal health
- Co-Biome's MetaXplore™ range provides three comprehensive gut microbiome and gastrointestinal health tests to support clinical decision making
- Using Co-Biome's MetaXplore™ range of tests combined with your clinical expertise can support informed prescribing and patient outcomes

CO-BIOME  
by MICROBA



50

## Introductory offer

- **30% discount** off the MetaXplore™ range available only to webinar attendees
- A discount link will be emailed to you
- The introductory offer will be available until 30 April 2023

### To gain referral access:

- Register: [co-biome.com/register/](https://co-biome.com/register/)




51

15 minutes

## Q&A from the chat

Hayley Parcell  
Dr Brad Leech  
Krystyna Sullivan



52

## The Co-Biome Team



**Hayley Parcell**  
Head of Co-Biome  
Healthcare  
*Nutritionist*



**Dr Brad Leech**  
Lead Clinical Educator  
*Integrated &  
Functional Medicine*



**Dr Paula Smith-Brown**  
Healthcare Science  
Liaison  
*Dietitian*



**Dr Alena Pribyl**  
Senior Scientist and  
Research Officer  
*Scientist*



**Krystyna Sullivan**  
Senior Clinical  
Application Specialist  
*Naturopath*



**Christal Jux**  
Clinical Application  
Specialist  
*Naturopath*



**Christine Stewart**  
Clinical Application  
Specialist  
*Nutritionist*



**Dr Angela Genoni**  
Project Scientist and  
Nutrition Lead  
*Nutritionist*



**Dr Nicola Angel**  
Head of Laboratory  
Operations  
*Scientist*



**Dr David Wood**  
Head of Bioinformatics  
Operations  
*Scientist*



**Dr Ken McGrath**  
Clinical Solutions  
Manager  
*Scientist*



**Alyssa Tait**  
Scientific Reviewer  
*Naturopath &  
Physiotherapist*



**Laima Hareer**  
Scientific Reviewer  
*Naturopath*



**Alana Krahe**  
Marketing Manager  
*Nutritionist*



**Jillian Foster**  
Content Specialist  
*Naturopath*



**Dr Areej Alsheikh**  
Bioinformatician  
*Scientist*

**CO-BIOME**  
by MICROBA

53

LIVE WEBINAR

# Thank you for attending

### Additional resources:

- MetaXplore™ Range Brochure
- How to refer handout
- Interpretation Guide (coming soon)
- Patient brochure (coming soon)

**CO-BIOME**  
by MICROBA

54

## References

- Ahola A.J., Lassenius M.J., Forsblom C., et al. Dietary patterns reflecting healthy food choices are associated with lower serum LPS activity. *Sci Rep* 7, 6511, 2017. <https://doi.org/10.1038/s41598-017-06885-7>
- Anhê F.F., Barra N.G., Cavallari J.F., et al. Metabolic endotoxemia is dictated by the type of lipopolysaccharide. *Cell Reports*, 36(11), 109691, 2021. <https://doi.org/10.1016/j.celrep.2021.109691>
- Asnicar F., Leeming E.R., Dimidi E., et al. Blue poo: Impact of gut transit time on the gut microbiome using a novel marker. *Gut*, 70(9), 1665–1674, 2021. <https://doi.org/10.1136/gutnl-2020-323877>
- Attaluri A., Jackson M., Valestin J., et al. Methanogenic flora is associated with altered colonic transit but not stool characteristics in constipation without IBS. *The American Journal of Gastroenterology*, 105(6), 1407–1411, 2010. <https://doi.org/10.1038/ajg.2009.655>
- Babidge W., Millard S., Roediger W. Sulfides impair short chain fatty acid beta-oxidation at acyl-CoA dehydrogenase level in colonocytes: Implications for ulcerative colitis. *Molecular and Cellular Biochemistry*, 181(1–2), 117–124, 1998. <https://doi.org/10.1023/a:1006838231432>
- Blachier F., Andriamihaja M., Larraufie P., et al. Production of hydrogen sulfide by the intestinal microbiota and epithelial cells and consequences for the colonic and rectal mucosa. *American Journal of Physiology. Gastrointestinal and Liver Physiology*, 320(2), G125–G135, 2021. <https://doi.org/10.1152/ajpgi.00261.2020>
- Byrd A.L., Liu M., Fujimura K.E., et al. Milieu Intérieur Consortium. Gut microbiome stability and dynamics in healthy donors and patients with non-gastrointestinal cancers. *J Exp Med*, 4:218(1):e20200606, Jan 2021. <https://doi.org/10.1084/jem.20200606>
- Chang Y., Deng Q., Zhang Z., et al. Glucagon-like peptide 2 attenuates intestinal mucosal barrier injury through the MLCK/pMLC signaling pathway in a piglet model. *Journal of Cellular Physiology*, 236(4), 3015–3032, 2021. <https://doi.org/10.1002/jcp.30068>
- de Mello V., Paananen J., Lindström J., et al. Indolepropionic acid and novel lipid metabolites are associated with a lower risk of type 2 diabetes in the Finnish Diabetes Prevention Study. *Sci Rep* 7, 46337, 2017. <https://doi.org/10.1038/srep46337>
- Derrien M., Alvarez A-S, de Vos W.M. The gut microbiota in the first decade of life. *Trends in Microbiology*, Vol 27, Issue 12, P997-1010, December 2019. <https://doi.org/10.1016/j.tim.2019.08.001>



## References

- Farhangi M.A., Vajdi M. Novel findings of the association between gut microbiota-derived metabolite trimethylamine n-oxide and inflammation: results from a systematic review and dose-response meta-analysis. *Critical Reviews in Food Science and Nutrition*, 60(16), 2801–2823, 2020. <https://doi.org/10.1080/10408398.2020.1770199>
- Gebrayel P., Nicco C., Al Khodor S., et al. Microbiota medicine: towards clinical revolution. *J Transl Med*, 7:20(1):111, Mar 2022. <https://doi.org/10.1186/s12967-022-03296-9>
- Gencer B., Li X.S., Gurmu Y., et al. Gut microbiota-dependent trimethylamine n-oxide and cardiovascular outcomes in patients with prior myocardial infarction: a nested case control study from the PEGASUS-TIMI 54 Trial. *Journal of the American Heart Association*, 9(10), e015331, 2020. <https://doi.org/10.1161/JAHA.119.015331>
- Ijssennagger N., Belzer C., Hooiveld G.J., et al. Gut microbiota facilitates dietary heme-induced epithelial hyperproliferation by opening the mucus barrier in colon. *Proceedings of the National Academy of Sciences of the United States of America*, 112(32), 10038–10043, 2015. <https://doi.org/10.1073/pnas.1507645112>
- Lyte J.M., Gabler N.K., Hollis J.H. Postprandial serum endotoxin in healthy humans is modulated by dietary fat in a randomized, controlled, cross-over study. *Lipids Health Dis* 15, 186, 2016. <https://doi.org/10.1186/s12944-016-0357-6>
- Matsuura M. Structural modifications of bacterial lipopolysaccharide that facilitate gram-negative bacteria evasion of host innate immunity. *Frontiers in Immunology*, 4, 109, 2013. <https://doi.org/10.3389/fimmu.2013.00109>
- McHardy I.H., Wu M., Shimizu-Cohen R., et al. Detection of intestinal protozoa in the clinical laboratory. *J. Clinical Microbiology*, Vol.52, No.3, 2014. <https://doi.org/10.1128/JCM.02877-13>
- Mokkala K., Houttu N., Koivuniemi E., et al. GlycA, a novel marker for low grade inflammation, reflects gut microbiome diversity and is more accurate than high sensitive CRP in reflecting metabolomic profile. *Metabolomics* 16, 76, 2020. <https://doi.org/10.1007/s11306-020-01695-x>
- Naveed A., Abdullah S. Impact of parasitic infection on human gut ecology and immune regulations. *transl med commun* 6, 11, 2021. <https://doi.org/10.1186/s41231-021-00091-4>
- Nighot M., Al-Sadi R., Guo S., et al. Lipopolysaccharide-Induced Increase in Intestinal Epithelial Tight Permeability Is Mediated by Toll-Like Receptor 4/Myeloid Differentiation Primary Response 88 (MyD88) Activation of Myosin Light Chain Kinase Expression. *The American Journal of Pathology*, 187(12), 2698–2710, 2017. <https://doi.org/10.1016/j.ajpath.2017.08.005>



## References

- Roager H.M, Hansen L.B.S, Bahl M.I, et al. Colonic transit time is related to bacterial metabolism and mucosal turnover in the gut. *Nature Microbiology*, 1(9), 16093, 2016. <https://doi.org/10.1038/nmicrobiol.2016.93>
- Schroemm A.B, Brandenburg K, Loppnow H, et al. Biological activities of lipopolysaccharides are determined by the shape of their lipid A portion. *European Journal of Biochemistry*, 267(7), 2000. <https://doi.org/10.1046/j.1432-1327.2000.01204.x>
- Sender R, Fuchs S, Milo R. Revised Estimates for the Number of Human and Bacteria Cells in the Body. *PLOS Biology*. 2016. <https://doi.org/10.1371/journal.pbio.1002533>
- Senthong V, Wang Z, Li X.S, et al. Intestinal microbiota-generated metabolite trimethylamine-n-oxide and 5-year mortality risk in stable coronary artery disease: the contributory role of intestinal microbiota in a COURAGE-like patient cohort. *Journal of the American Heart Association*, 5(6), 2016. <https://doi.org/10.1161/JAHA.115.002816>
- Tuomainen M, Lindström J, Lehtonen M, et al. Associations of serum indolepropionic acid, a gut microbiota metabolite, with type 2 diabetes and low-grade inflammation in high-risk individuals. *Nutr Diabetes*. 25;8(1):35, May 2018.. <https://doi.org/10.1038/s41387-018-0046-9>
- Whipps J.M, Lewis K, Cooke R.C. Mycoparasitism and plant disease control 161–187. In: Burge, NM (editor), *Fungi in Biological Control Systems*. Manchester University Press; 1988. P. 176.
- Yang J, Guo Y, Lee R, et al. Pomegranate Metabolites Impact Tryptophan Metabolism in Humans and Mice. *Curr Dev Nutr*. 25;4(11):nzaa165, Nov 2020.. <https://doi.org/10.1093/cdn/nzaa165>
- Zamyatina A, Heine H. Lipopolysaccharide Recognition in the Crossroads of TLR4 and Caspase-4/11 Mediated Inflammatory Pathways. *Frontiers in Immunology*, 11, 585146. 2020. <https://doi.org/10.3389/fimmu.2020.585146>
- Zhernakova A, Kurilshikov A, Bonder MJ, et al. Population-based metagenomics analysis reveals markers for gut microbiome composition and diversity. *Science*. Apr 29;352(6285):565-9, 2016. <https://doi.org/10.1126/science.aad3369>